

WWW.NORABODE.COM



NORA BODE

# OXYJET GO

NATURALLY BEAUTIFUL SKIN WITHOUT NEEDLES:  
GENTLE, EFFECTIVE AND SAFE!

THE GLOW! ONE.TWO.GO . GET THE GLOW! ONE.T





NORA BODE OXYJET

*The concept*

## MODE OF ACTION

*The OXYJET System is a high tech oxygen treatment for your beauty needs. For over 20 years, NORA BODE Kosmetik have specialised in the needle free oxygen pressure injection and are the original inventors and patent holders of this method.*

*In cosmetic and medical practices worldwide, the Bode OXYJET is the method of choice when it comes to gentle but effective treatment of skin problems, especially wrinkles, pigment irregularities and loss of elasticity.*

NORA BODE now introduces another unique innovation into the World of Beauty, which means you can supply your skin with concentrated active ingredients and molecular encapsulated oxygen at home or wherever you please; smoothing out wrinkles and revitalizing your skin – gentle, effective and safe!

Traditionally, "Made in Germany" and NORA BODE Kosmetik manufacturing assures all of our customers of the outstanding quality of all our products. The new OXYJET GO system was developed and manufactured in Germany – producing you a high quality and exclusive product.

**NORA BODE KOSMETIK**




# INTRODUCTION

***This is how the OXYJET treatment works***

## ***The professional treatment in the practice***

*Using pulsed oxygen pressure, special cosmetic formulations  
and pure oxygen are "shot" into the deep layers of the skin -  
without needles - gentle, effective and safe!*




*The professional treatment in the practice.*

*Oxygen massage against swellings.*



*Pulsed oxygen pressure injection.*




The intensity of the pulsed pressure is fixed and this enables the active ingredients to penetrate through the natural interstices between the epidermal skin cells into the deep layers of the skin and build up a deposit there. An additional energy boost is provided by the concentrated oxygen used in the treatment as a carrier for the active ingredients. The Fraunhofer Institute for Biomedical Engineering in St. Ingbert, Germany, has impressively confirmed the strong drug-penetration-accelerating effect of the Bode OXYJET in test series. Different from needle injections, all cosmetic problem areas can be successfully treated with the OXYJET: The whole face, especially the delicate skin around the eyes and mouth, neck, décolleté and body.


5-10 professional treatments are necessary to achieve optimal and long lasting results. To keep the deposit of active ingredients in the skin as long as possible and to stabilize the treatment result, it is essential to use the especially for the OXYJET treatment designed NORA BODE aftercare products. The result of a professional OXYJET treatment course: Smooth, more regular pigmented, fresh skin with significantly reduced wrinkle depth and a visible glow.

# THE OXYJET GO TREATMENT

*in the comfort of your home*



Before the treatment



After the treatment

## High pressure for beautiful skin

*Supply your skin with microencapsulated oxygen and concentrated active ingredients, smooth out lines, reduce pigmentation irregularities and give your skin a fresh look, whenever and wherever you like: This can all be achieved with the new OXYJET GO system by NORA BODE.*

Modified for exclusive use at home or on the move, you are now able to treat the problem areas of your skin systematically. You apply the special preparations included in the set and then push them gently into the deep layers of the skin using pulsed pressure. In addition, you can perform a treatment against swellings around the eye area.



Dermatest, Münster, Germany has tested the OXYJET GO regarding its efficiency. The results after only one treatment:

- The OXYJET GO pulses active ingredients deep into the skin- without needles!
- The epidermal intercellular spaces are filled up, the volume of the skin increases.
- A deposit of active ingredients is formed.
- The skin remains intact and unhurt.

## What happens to your skin


You apply the OXYJET 1 and 2 preparations especially developed for the treatment and speed up their penetration using the OXYJET GO pulsed pressure device. In this way you can build up a deposit of active ingredients in your skin, which provides an immediate smoothing effect as well as a long term regenerating effect on your skin appearance.

**What you can expect**

The key ingredient "molecular oxygen" can help your skin to look fresher, smoother and radiant at any age.

Already in the early 2000s scientists of the Ruhr University Bochum have discovered that the skin is exclusively supplied by oxygen from the outside- down to the dermal layers! In addition, studies show that increased supply of molecular oxygen stimulates collagen formation and replenishes the energy stores of the skin.

The NORA BODE OXYJET 1 serum contains an active ingredient also used in medicine, which releases molecular oxygen to the skin, so that by the application of the pulsed pressure treatment it can penetrate deep into your skin and unfold its effect there.



**The OXYJET GO device**

Equipped with two treatment functions and modern monitoring technology, the device offers maximum treatment comfort:

**PULSE-MODUS**

By precisely defined pressure surges, the special preparations OXYJET 1 and 2 are pulsed deep into the skin.

**FLOW-MODUS**

By application of the constant flow, you can gently massage your skin and soften puffiness around the eye area.

**LIGHT RING INDICATOR**

The light ring in the device head indicates the current operating status of your OXYJET GO device. In SELECT mode and during treatment the light ring illuminates constantly blue. After the maximum treatment time the ring illuminates constantly red.

**MULTIPLE PROTECTION SYSTEM**

The integrated multi-level safety system MPS (Multiple Protection System) ensures maximum device safety and supports the user with optical and acoustic signals.



The ring illuminates blue and red.



**The NORA BODE treatment preparations**

*You will find two special preparations for use with the OXYJET GO system in your set.*

*They contain the necessary microencapsulated oxygen and further high-tech active ingredients*

*essential for a successful treatment and a radiant skin appearance. Both are free*

*from parabens, fragrances, alcohol and PEGs and essential for the application with the*

*pulsed pressure.*

**OXYJET 1 SERUM**

This special preparation releases molecular oxygen and is the basis of your treatment. The active ingredient used here can demonstrably reduce wrinkles and smooth the skin. In addition, pigment-regulating and intensively moisturizing agents are included.

**Perfluorodecaline** is a substance similar to the red blood pigment hemoglobin, able to bind oxygen and to release it at the same time. By its encapsulation, Perfluorodecaline penetrates deep into the skin and there releases the oxygen to the natural intercellular spaces, so that the skin cells can absorb it easily. So the skin can be supplied with energy, the formation of new cells be increased and by that the skin be vitalized and smoothed.

**ATP** stores the energy gained from the oxygen so that it is available to the skin cells for a longer period of time.

**Hyaluronic Acid** in form of different molecular weights penetrates deep into the skin, forms a deposit there and plumps up the skin from the inside out. Throughout the whole day the skin is supplied with new elasticity, firmness and tone.

**Meiritage®** consists of root extracts of various Asian herbs. The active substance-complex intensively moisturizes, stimulates the collagen synthesis and can provide an evenly pigmented complexion.

**Squalane and Glycosphingolipids** fill up the skin's own lipids, regenerate the cell membranes and provide a soft and supple skin feeling.

**OXYJET 2 SERUM**

This special preparation aims to give your skin new volume from the inside out. A special plant ingredient approved in Ayurvedic Medicine over hundreds of years can enhance the plumpness of the subcutis naturally. Additionally tightening, soothing and eye zone deswelling ingredients are included.

**Commiphora Mukul Extract** stimulates the synthesis of lipids deep in the skin and by that can enhance the natural plumpness which usually only exists in young skin. Applied with the OXYJET GO the active ingredients can penetrate deep into the skin and gradually unfold its filling effect – fully naturally and without needles.

**Sweet Almond Oil** contains a high content of unsaturated fatty acids and vitamins and therefore nurtures the skin without being comedonic.

**Glycine Soy Oil** is rich in unsaturated omega-6 fatty acids, strengthens the barrier function of the skin and provides a velvety soft skin feel.

**Pullulan** is a natural saccharide which forms a tightening film on the face and by that provides clearer face contours.

**Acetyl Tetrapeptide-5** is a synthetically derived active ingredient that can reduce puffiness and dark discoloration around the eyes and provide a more youthful appearance especially in this area.

## THE TREATMENT

### „step by step“

*Your OXYJET GO device is designed for the treatment of your complete face. Split your face into areas for best results: Forehead, eye and mouth area are treated with pulsed pressure step by step.*

As we age we all suffer from tired-looking skin with wrinkles around the eyes and mouth usually among the first signs. Furthermore, pigment irregularities and decreasing elasticity occur at a young age with this just being a sign of normal living. However, you can do something quite naturally to enable you to look fresher, smoother and younger. The OXYJET pulse pressure treatment is gentle, effective and safe for your naturally fresh appearance.

### **The preparation of your skin**


As usual, cleanse your face with a mild cleansing gel, foam or milk. Dry your skin carefully afterwards.





**Step by step - the treatment**

Take 4 to 5 drops of the OXYJET 1 serum, spread it thoroughly over your complete face and pad in gently. This way you can supply your skin with oxygen, energy-giving and pigment-regulating active ingredients, thus creating the basis for every OXYJET GO treatment.



*Distribution of OXYJET 1 serum.*



*Distribution of OXYJET 2 serum.*



*Treatment of the eye area.*



*Treatment of the eye area.*

Then start with the treatment of the face area you have chosen, here for example the lower eye zone: Spread 1-2 drops of OXYJET 2 serum below the eye and pulse it deep into the skin by the application of the OXYJET GO device.

If you suffer from puffiness in the lower eye area, you can perform a gentle, decongestive massage using the "Flow" mode. For the treatment of forehead and mouth area please proceed in the same way.

### **How often you should use the OXYJET GO device**

---

In order to build up a deposit of active ingredients in your skin and to get a fresher, smoother and more even complexion, you should use your device 3 times a week preferably in the evening for a period of 4 weeks. There should always be at least 1 treatment-free day between the single treatments. When you have finished the 4-week treatment cycle, you should use your OXYJET GO device once a week in order to maintain the deposit of active ingredients in your skin. For an immediate effect, if you want to freshen up quickly for a special occasion, you can use the device at any time in between.

## AFTER THE TREATMENT

---

### ***The final care of your skin***


After you have finished the treatment with your OXYJET GO device, you should provide your skin with your daily day/night care. Both NORA BODE treatment supplements are highly concentrated active substance serums that you should use during your treatment cycle every day in the evening under your usual care. In this way you support the deposit of active ingredients created by the OXYJET GO treatment and thus ensure an optimal treatment result.



*Since more than 25 years medical and cosmetic practises worldwide use NORA BODE Kosmetik,  
our OXYJET devices and preparations successfully.*

## LEARN MORE ONLINE

---



**For more information on NORA BODE Kosmetik  
and the pulsed oxygen pressure injection, visit  
[www.norabode.com](http://www.norabode.com).**

**OXYJET treatment preparations**

Please contact the authorized OXYJET GO distributor in  
your country.

**Spare parts like nozzles and power adapter**

Please contact the authorized OXYJET GO distributor in  
your country.

[www.norabode.com](http://www.norabode.com)





THE GLOW! ONE·TWO·GO · GET THE GLOW! ONE·T

